



C2000-1/2/4
Semi-Auto Coagulation Analyzer

Dual magnetic circuit beads method
— A promise of reliable result

C2000-1/2/4 Semi-Auto Coagulation Analyzer

Mechanical detection method, insensitive to abnormal samples that are hemolysis, lipemic, icteric etc.



Technical specifications

Testing method	Coagulation method (mechanical)
Testing items	PT, APTT, TT, FIB etc.
Sample volume	≤200μL
Reagent volume	50μL-200μL
Timing start	Electronic Start-Pipette
Channel consistency	4%

Accuracy	±3%
Testing Channel	1/2/4
Incubation Channel	4/6/16
Reagent Position	2/3/4
Printing	Built-in Printer for C2000-2/4
Optional	Mini printer for C2000-1